## VITA VIONIC® WAX fräsbarer Rohling aus Synthetikwachs

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### Instructions for use

## VITA VIONIC® WAX milled synthetic wax blank

Please read the instructions for use before using the product.

### Intended purpose

VITA VIONIC WAX is a blank made of millable and dimensionally stable synthetic wax with a high melting temperature for the fabrication of full-sized wax try-ins and wax-ups.

### Patient target group

No restrictions

### Intended user

Dental professionals only: dentists and dental technicians

### Indication

Indication range:

- VITA VIONIC WAX White is used for the fabrication of full-sized try-ins. They are used to test the midline, occlusion plane and phonetics.
- VITA VIONIC WAX Pink is used for the fabrication of denture bases for try-in. For this purpose, denture teeth are fixed in the milled cavities. All parameters can be checked, and if necessary, corrections can be implemented.

### Contraindication

- Not suitable for the fabrication of definitive dental restorations.

### Product application/control

#### Workflow:

VITA VIONIC WAX White:

1. Clamp the wax disc.
2. Mill the wax disc (maxilla and mandible),
3. Separate the try-in,
4. Check the try-in against the master model,
5. Intraoral try-in,
6. Implement any required changes digitally using the CAD software,
7. If necessary, perform another try-in.

VITA VIONIC WAX Pink:

1. Clamp the wax disc.
2. Mill the wax disc (maxilla and mandible),
3. Separate the try-in,
4. Use wax to fix denture teeth (for conventional denture teeth, prior CAM modification is required),
5. Check against the master model,
6. Intraoral try-in,
7. Scan in any required changes as a new registration, and implement them digitally using the CAD software,
8. Alternative conventional method of finalizing the dental restoration.

The product must be processed at room temperature. Wet machining is recommended for the milling process.

#### System requirements:

For the VITA VIONIC WAX system requirements, see [www.vita-zahnfabrik.com/vionicwax](http://www.vita-zahnfabrik.com/vionicwax).

Detailed information on processing can be found in the document "Instructions for Use I full version" at [www.vita-zahnfabrik.com/vionicwax](http://www.vita-zahnfabrik.com/vionicwax).

### Recommended product combinations

For recommended product combinations with the VITA system components and corresponding processing recommendations and restrictions regarding the combination possibilities of the product, please refer to the Instructions for Use | full version.

### Storage/disposal

Store at room temperature. Protect from direct sunlight and temperatures above 25°C.

The countersigned products labelled with a pictogram for hazardous substances are to be disposed of as hazardous waste. Recyclable waste (such as attachments, paper, plastics) must be disposed of using appropriate recycling systems. If necessary, contaminated product residues should be pretreated in accordance with regional regulations and disposed of separately.

### Technical data

Physical properties	Unit of measure	Value
Solidification point	°C	83
Drop point	°C	89
Flash point	°C	>140
Penetration at 25 °C	0.1 mm	9
Viscosity at 120 °C	mPas	9.7 - 10.1

Chemical composition	Wt%
Paraffin waxes	65 - 85
Pentaerythrityl	10 - 30
Polyethylene	5 - 12

The technical data are averages and may fluctuate slightly between delivered batches.

### Product reliability

Information on reporting serious incidents in connection with medical devices, general risks associated with dental treatments, residual risks and (if applicable) short clinical safety and performance reports (SSCPs) can be found at [https://www.vita-zahnfabrik.com/product\\_safety](https://www.vita-zahnfabrik.com/product_safety).

The safety data sheets can be downloaded at [www.vita-zahnfabrik.com/vionicwax](http://www.vita-zahnfabrik.com/vionicwax) or requested by fax at (+49) 7761-562-233.

### Disclaimer

Please note: Our products must be used in accordance with the instructions for use. We do not accept liability for any damage resulting from incorrect handling or usage. The user is furthermore obliged to check the product before use with regard to its suitability for the intended area of applications. We cannot accept any liability if the product is used in conjunction with materials and equipment from other manufacturers that are not compatible or not authorized for use with our product and this results in damage. The VITA Modulbox is not necessarily a component of the product. Date of issue of these instructions for use: 2022-03

After the publication of this information for use, any previous versions become obsolete. The current version can be found at [www.vita-zahnfabrik.com](http://www.vita-zahnfabrik.com).

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Instructions for use | short version

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VITA Zahnfabrik has been certified, and the following products bear the mark:

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VITA VIONIC® WAX

UK UK Responsible Person: Obelis UK Ltd.  
CA Sandford Gate, Oxford, OX4 6LB, United Kingdom

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CH REP VITA Zahnfabrik H. Rauter GmbH & Co.KG, Bad Säckingen (Germany)
Zweigniederlassung Basel c/o Perrig AG, Max Kämpf-Platz 1, 4058 Basel
VITA Zahnfabrik H. Rauter GmbH & Co.KG
Gebdi Dentalproducts GmbH