# HOJA DE INSTRUCCIONES

es

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# INSTRUCTION SHEET

en

G-CAM. Graphene nanoreinforced biopolymer disc for CAD/CAM drilling. Dental use.

Please read this instruction sheet carefully before using the product.

## 1. What is G-CAM?

G-CAM is a thermostable acrylic with a principal base of polymethyl methacrylate resin doped with graphene (allotropic form of carbon), suitable for the creation of dental prostheses using CAD/CAM technology.

The product is available in 2 different formats:

- **G-CAM MONOCHROMA.** Monochrome colour disc format for full anatomical monolithic restorations or aesthetic coated meso-structures.
- **G-CAM MULTICHROMA.** Multichromatic disc format for performing full anatomical monolithic restorations.

The disc dimensions are as follows: 98.5 mm (universal anchorage) or 95 mm (Zirkonzhan anchorage) in diameter; 14, 16, 18, 20, 22, 24, 26, 28 and 30 mm.

## 2. Indications

G-CAM has been designed to provide an integral solution for dental prostheses using the CAD/CAM technique.

## 3. Contraindications

Although the elimination of residual monomers which can cause allergies is tested during the manufacturing method, there are no recognised contraindications, except in the cases of adverse reactions to polymethyl methacrylate.

## 4. Risk and safety indications

### 4.1 Working with CAM

The G-CAM disc can be used either wet or dry. In case of dry operation, it will be necessary to suction any debris resulting from the drilling operation. The drills may break if improper tools are used. Although graphene helps to dissipate the heating of materials, an inappropriate strategy may cause the overheating and deformation of the material.

### 4.2 Work in the laboratory

The occupational safety measures indicated in the regulations which are currently in force will be applicable. A protective mask and suction equipment must be used to prevent the inhalation of dust. Safety goggles should be worn to avoid damage to the eyes caused by foreign particles. The information and advice contained in this instruction sheet should be adhered to in order to safeguard the quality and safety of the material.

### 5. Indications regarding storage and expiration

The product should be stored at room temperature in its original packaging.

### 6. Processing and application

The minimum thickness of the walls of the reconstructions should be 0.6 mm in occlusal and 0.4 mm in cervical areas. The cross-sections of connectors in the anterior area should not be less than 9 mm<sup>2</sup> or less than 13 mm<sup>2</sup> in posterior areas. There should be no more than three elements between the two pillar pieces. Review all the G-CAM design parameters established for all the different dental treatments.

G-CAM can only be used in drilling systems without RFID connection (radio frequency identification), as our product is designed for completely open systems.

The PMMA drilling strategy is used.

### 7. Technical characteristics

- Elastic modulus: 3200 ± 7% MPa
- Bending strength: 140 ± 7% MPa
- Surface hardness: 88 Shore
- Water absorption: 4 µg/mm<sup>2</sup>
- Residual monomer: <0,004 %

### 8. Other observations

G-CAM must be used only by dentists, dental technicians and qualified personnel with specific training. Only original tools and parts may be used for processing.

This instruction sheet corresponds to the current technique and our own experience. The product may only be applied for the indications described in point 2. The user is solely responsible for the application of the product. As manufacturers, we do not assume any responsibility for failed results, given that we do not participate in the process. Any claims for damages are limited exclusively to the commercial value of our products.

### 9. Applied symbols

|  |  |
|--|--|
|  | Lot number   |
|  | Reference number   |
|  | Colour   |
|  | Protect from sunlight  |
|  | Check the instructions of use                                    |
|  | Manufacturer   |
|  | Take into account the documentation that accompanies the product |
|  | Format   |
|  | Store in a dry area  |
|  | Diameter   |
|  | Thickness  |
|  | Do not reuse   |
|  | 1984<br>CE: EC Certification<br>1984: Notified Organism KIWA     |

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# Design parameters

## G-CAM disc

|                      | Thicknesses   | Occlusion   | Pontic  | Connector cross section   | Submerged implant   | Vertical dimension   | Milled strategy |
|----------------------|---|---|---|---|---|--|-----------------|
| <b>Crown</b>         | Cervical 0,4 mm<br>Occlusal 0,6MM   | Light occlusion in intercuspal occlusion and without contacts in excursive movements.   |   |   |   |  | PMMA            |
| <b>Bridge</b>        |   | Light occlusion in intercuspal occlusion and without contacts in excursive movements.   | 3 pontics between crowns or abutments.  | Anterior: 9 mm <sup>2</sup> or more.<br>Posterior: 13 mm <sup>2</sup> or more                       |   |  | PMMA            |
| <b>Meriland</b>      | Minimum 0,6 mm  | No occlusion, recommended for <b>PROVISIONAL TREATMENTS</b> .   | 1 pontic with supports on both sides. Only in anterior area.                                  | Minimum 9 mm <sup>2</sup>   |   |  | PMMA            |
| <b>Inlays onlays</b> | "Cervical 0,4 mm<br>Occlusal 0,6 mm"  | Light occlusion in intercuspal occlusion and without contacts in excursive movements.   |   |   |   |  | PMMA            |
| <b>Veneers</b>       | Minimum 0,3 mm  | Upper veneers singulum natural occlusion with the lower arch and lower veneers soft occlusion and no interference in protrusive and lateral movements.  |   |   |   |  | PMMA            |
| <b>Implants</b>      | Radius abutment base: 0,15 mm and the thickness around the abutment should be approximately 2,5 mm. | <b>Crowns:</b> light occlusion in intercuspal occlusion and without contacts in excursive movements.<br><b>Half arch:</b> functional occlusion in intercuspal occlusion and bilateral group function.<br><b>Full arch:</b> balanced articulation. | 3 Pontic between abutment. If we have more than 3 pontics, make assessment of load factors. * | Anterior: 9 mm <sup>2</sup> or more.<br>Posterior: 13 mm <sup>2</sup> or more                       | Place with an appropriate height abutment, to raise the connection up to the gingiva as possible and thus ensure the thickness of the emergency profile in the restoration and avoid descentration. | If the height of the tibase do not exceed the 50% of the height of the restoration, it will have to be replaced with abutment, to ensure the stability of the structure and avoid descentration. | PMMA            |
| <b>Cantilever</b>    |   | Without occlusion   | 1 pontic, maximum 11 mm in length. *  | Minimum 13 mm <sup>2</sup> and do not individualize the structure at the lingual or palatinal area. |   |  | PMMA            |

\* If there are more than three pontics between abutments it will be necessary to evaluate the load factors, but in case of doubt we recommend reinforce internally the G-CAM with a structural material.

\* If there are more than 1 cantilever bridge, it is necessary to reinforce internally the G-CAM with a structural material.

Made in Spain

|                                      |                                |                          |  |
|--------------------------------------|--------------------------------|--------------------------|--|
| <b>Crowns</b>                        | <b>Veneers</b>                 | <b>Submerged Implant</b> | <b>Vertical dimension</b>  |
|                                      |                                |                          |  |
| <b>Thickness around the abutment</b> | <b>Connector cross section</b> | <b>Up to 3 pontics</b>   | <b>Cantilever</b>  |
|                                      |                                |                          |  |
| <b>Onlays / Inlays</b>               | <b>Milling strategy</b>        | <b>Finishing</b>         |  |
|                                      |                                |                          | <p><b>G-CAM</b></p> <p>Graphene nanoreinforced biopolymer disc for CAD/CAM drilling</p> <p>Made in Spain</p> |